

LUCAS™2 BATTERY CHARGER

Instructions for Use

GB

LUCAS™2 Battery Charger (also referred to as the Charger) is an accessory to the LUCAS™2 Chest Compression System. It is a stand-alone unit intended to charge the LUCAS™2 Battery (also referred to as the Battery).

DECLARATION OF CONFORMITY




LUCAS™2 Battery Charger complies with the requirements of the following European Directives:



2006/95/EC Low Voltage Directive (LVD)
2004/108/EC EMC Directive

The device is marked with the CE-symbol: **CE 0434**

© Copyright JOLIFE AB 2009. All rights reserved.
LUCAS™ is a trademark of JOLIFE AB.

READ THESE INSTRUCTIONS BEFORE USING THE LUCAS™2 BATTERY CHARGER:

	Caution – see instructions for use All operators must read the complete Instructions for Use before using LUCAS™2 Battery Charger.
	LUCAS™2 Battery Charger is only to be used with JOLIFE products.
	The mains socket should be easy accessible. If an operational error occurs, the plug should be immediately removed from the socket.
	LUCAS™2 Battery Charger is made for indoor use and should not be exposed to water or dust.
	LUCAS™2 Battery Charger contains dangerous voltage and the cover should not be removed. All service or maintenance work should be carried out by a qualified technician who can get assistance by contacting the manufacturer's agent.
	WARNING The LUCAS™2 Battery Charger has internal fuses which blow if a fault occurs in the charger. Such faults must be repaired by a qualified service technician.
	Check that the mains cord is not damaged. If the mains cord is damaged, the LUCAS™2 Battery Charger must not be used.
	Do not cover up LUCAS™2 Battery Charger when in use.
	The LUCAS™2 Battery Charger is turned on by connecting it to the mains socket. Disconnecting it from the mains socket turns it off.

	Electrical and electronic equipment must not be discarded with domestic waste, but be disposed of using separate collection, treatment, recovery/recycling and environmentally safe methods.
	Year of manufacture is printed on the type label next to this symbol.





PREPARATION FOR USE

Connect the mains cord to the outlet on the Charger.

CHARGING INSTRUCTIONS

1. Place a LUCAS™2 Battery in the Charger.
2. Connect the mains cord to the mains socket. The LED (Light-Emitting Diode) on the top of the Charger shows the Charger status:
 - a. Red – Charging
 - b. Yellow – Charging (80-95% completed)
 - c. Green – Fully charged / standby
3. When the Battery is fully charged (LED is green):
 - a. Disconnect the mains cord from the mains socket
 - b. Remove the Battery from the Charger

TECHNICAL SPECIFICATIONS:

Type	Lithium-Ion Battery Charger		
Number of cells	7		
Input	90-264 VAC, 47-63 Hz, max 0.9A		
Temperature range – operating	-20 °C – +40 °C		
Temperature range – storage	-25 °C – +85 °C		
Insulation class	II		
Red	Charging	LUCAS™2 Battery charge level: 	1.2 A , 29.4 V
Yellow	Charging (80-95% completed)	LUCAS™2 Battery charge level: 	29.4 V , < 0.5 A
Green	Ready / standby	LUCAS™2 Battery charge level: 	

For more information - please contact the manufacturer Jolife AB

JOLIFE AB, Scheelevägen 17, SE-223 70 Lund, Sweden Tel: +46 (0) 46 286 50 00, Fax: +46 (0) 46 286 50 10, Email: info@jolife.com
jolife.com

LUCAS is a trademark of JOLIFE AB.
CE 0434 © 2009 JOLIFE AB