

LUCAS[®] 2 CHEST COMPRESSION SYSTEM

LUCAS 2 and Larger Patients



Weight: 180 lb
Height: 6 ft



Weight: 330 lb
Height: 6 ft, 5 in



Weight: 320 lb
Height: 5 ft, 10 in





Proper Placement is on Upper Chest

Many large patients carry their weight through their mid-section. When LUCAS is properly applied, it rests over the upper chest or sternum, with enough room for the compressor to give 2 inch compressions. As long as the unit snaps into the backboard and the suction cup does not compress the patient's chest in the start position, it will operate as intended.

LUCAS 2 CHEST COMPRESSION SYSTEM



Patients Eligible for Treatment:

- Sternum height of 6.7 to 11.9 inches (17 to 30.3 cm)
- Maximum chest width: 17.7 inches (45 cm)

The use of LUCAS is not restricted by patient weight.

For further information please contact your local Physio-Control representative or visit our website at www.physio-control.com



**Physio-Control
Headquarters**
11811 Willows Road NE
Redmond, WA 98052 USA
Tel 425 867 4000
Fax 425 867 4121
www.physio-control.com

**Physio-Control
Operations
Netherlands B.V.**
Keizersgracht 125-127,
1015CJ Amsterdam NL
Tel +31 (0)20 7070560
Fax +31 (0)20 3301194
www.physio-control.nl

**Physio-Control
UK Sales Ltd**
Concorde House
Trinity Park, Solihull
Birmingham
B37 7UQ
United Kingdom
Tel/Fax 0808 258 0094 option 1

**Physio-Control
Australia Pty Ltd**
Suite 4.01
15 Orion Road
Lane Cove
NSW 2066
Australia
Toll Free Tel 1800 987 982
Toll Free Fax 1800 890 892



Jolife AB, Scheelevägen 17, SE-223 70 LUND, Sweden